
#### **Aseptic Technique**

- Refers to use of methods & procedures that prevent microbial contamination while performing a given task
- Refers to procedures performed in:
  - Hospital (operating theatres, nursing care)
  - o Research laboratories (e.g. working with cell culture)
  - Pharmacy: Industry (in production of parenteral dosage forms, ophthalmic preparation etc.)
     & Hospital pharmacy (preparation of parental solutions)
- Requires use of personal protective equipment (PPE): Disposable latex/nitrile gloves, Lab coat,
   Safety glasses, Face mask, Hair net & Appropriate footwear
- Use of PPE is essential to protect both personnel & product

#### **Aseptic Workstations**

## **Laminar Flow Hood** Air flows in one direction without turbulence HEPA (High Efficiency Particulate Air) filter: Removes particles >0.3 μm Provide protection only to product No protection for personnel Class I Biosafety Open fronted Cabinets Negative pressure systems Provide personal protection & environmental protection With HEPA filtration system fitted Do not provide product protection Class II Biosafety Negative pressure system ventilated cabinet: Provides personnel **Cabinets** protection HEPA-filtered air circulated into cabinet: Provides product protection Exhaust air is HEPA filtered: Provides environmental protection



### **Disinfectants & Antiseptics**

### **Characteristics of Disinfectants, Antiseptics & Preservatives**

	Disinfectant	Antiseptic	Preservative	
	Used on surface of inanimate objects (e.g. table, floor)	Used on living tissue     (e.g. skin, mucous     membranes)	<ul> <li>Not disinfectant/ antiseptic</li> <li>Additives that prevent microbial growth (bacteriostatic)</li> </ul>	
Activity	Biocidal	Biocidal	Biostatic	
Toxicity	Can be toxic	Non-Toxic	Non-Toxic	
Irritant	Can be an irritant	Non-Irritant	Non-Irritant	
Stability	Can be kept in stable form until required	Can be kept in stable form until required	Needs a long shelf-life	
Reactivity	Non-reactive with non-living items	Non-reactive with live tissues	Non-reactive with other ingredients	

#### **Levels of Disinfection**

	Low	High
Microorganisms	Most vegetative bacteria	All microorganisms
Killed	Most viruses	Most bacterial spores
	Most fungi	
Microorganisms	Mycobacteria	Some bacterial spores
Surviving	Bacterial spores	• Prions
	• Prions	
When to Use	Non-critical medical devices (e.g.	Critical medical devices
	wheelchair, blood pressure cuffs)	

# **Disinfectant Types**

Hospital	Bactericidal efficacy is mandatory (excluding Mycobacteria)		
Disinfectants			
Disiniectants	Must be general purpose suitable for many surfaces		
	Alcohols & phenols are not high-level disinfectants unlike hydrogen peroxide,		
	glutaraldehyde		
	<ul> <li>Hospital grade disinfectant should be able to kill common hospital pathogens</li> </ul>		
	MICROBIOLOGICALLY TESTED:		
	has been <i>proven</i> to kill:		
	✓ Staphylococcus aureus (MRSA ✓ Salmonella choleraesuis		
	or Golden Staph) ✓ Proteus vulgaris		
	✓ Pseudomonas aeruginosa ✓ Influenza virus		
	✓ E coli ✓ Herpes Simplex virus		
	✓ Enterococcus Faecalis (VRE) ✓ Hepatitis B Group virus ✓ Acinetobacter		
	✓ Acetobacter		
	- Acetobactes		
Household	No mandatory levels		
Disinfectants	General purpose		

#### **Choosing Right Antimicrobial Agent**

Toxicity	<ul> <li>Particularly important for antiseptic &amp; preservatives</li> </ul>	
	When is chemical toxic? (MSDS)	
	Everything is toxic at high enough conc.	
	Should be considered for disinfectants:	
	<ul><li>Is it volatile? Will inhalation be an issue?</li></ul>	
	<ul> <li>What personal protective equipment will be required?</li> </ul>	
Properties of	• Is it a solvent? Will it stain/corrode? Are these acceptable?	
<b>Chemical Agent</b>	What is rate & extent of microbial killing over time?	
	Is it affected by temperature/pH?	
	• Does it have a negative effect on living tissue? If so, it cannot be an	
	antiseptic/preservative.	
	Is it toxic once ingested? If so, it cannot be a preservative	
Intended	How long does it require to be effective?	
Application	Is it compatible with intended surface?	
	• Is it compatible with other formulation ingredients? Especially important	
	for preservatives. (e.g. not using bleach with hot water)	
	Is it still effective at non-toxic conc?	
	Are there any long-term effects?	

#### Microbial What is the bioburden? (Bioburden is level of microbial contamination) Challenge What types of organisms are present? Is there a pyrogen present? Cysts Prion Coccidia Bacterial spores Mycobacteria Small enveloped Lipid enveloped Gram-positive Large enveloped Gram-negative Fungi viruses bacteria bacteria **Prions** Transmissible infective proteins Exist in two forms (native form PrP<sup>C</sup> & misfolded infective form PrPSc) PrPSc can refold native/healthy PrPC Extremely stable **Bacterial** Gram+ bacteria (e.g. Clostridium tetani & C. botulinum) **Endospores** form spores Very resistant to environmental factors Bacteria can form endospores in 6-8 hours Remain viable for a long period of time Form a vegetative cell in 90 mins in favourable conditions Asymmetric division Germination Spore release Engulfment Coat formation Cortex synthesis **Environmental** Item must be cleaned before it can be disinfected **Factors** Organic matter disrupts antimicrobial effect as organic matter protects microbes, adsorbs disinfectant & chemically inactivates disinfectant Bodily fluids can prevent affect disinfection Biofilms can tolerate disinfectants